



PRESCRIPTION

PATIENT ACCESS

PREPARATION

APPLICATION

DISCHARGE
MANAGEMENT

Original Perfusor® Line

THE CHOICE IS YOURS

Original Perfusor® Line

Color coded tubes for an individual assignment



POLYETHYLENE
(PE)



NOT MANUFACTURED
WITH DEHP



NOT MANUFACTURED
WITH PVC



2 BAR PRESSURE
RESISTANT

- Option to assign high potential drug by color
- Quick identification of specific drug applications
- Designed to reduce risk of confusion
- Low priming volume
- Excellent infusion start-up characteristics in combination with syringe pumps
- 3 different color coded tubes in blue, magenta, green

| Product | Color | Material | Connection | Tube length (cm) | Diameter (mm) | Priming volume (ml) | Units per box | Code no. (REF) |
|-------------------------|---------|---------------|------------|------------------|---------------|---------------------|---------------|----------------|
| Original Perfusor® Line | Blue | PE (PVC-free) | Luer-Lock | 150 | 1.0 x 2.0 | 1.27 | 100 | 8723020 |
| | | | | 200 | | 1.67 | | 8723021 |
| | Magenta | PE (PVC-free) | Luer-Lock | 150 | 1.0 x 2.0 | 1.27 | 100 | 8723025 |
| | | | | 200 | | 1.67 | | 8723026 |
| | Green | PE (PVC-free) | Luer-Lock | 150 | 1.0 x 2.0 | 1.27 | 100 | 8723023 |
| | | | | 200 | | 1.67 | | 8723024 |