

# Easypump® II – Product Specifications

## LONG-TERM INFUSION

Designation	Color Code	Code No.* (REF)
Easypump® II LT 60-12-S	5 ml/h	4540002
Easypump® II LT 500-12.5-S	40 ml/h	4540003
Easypump® II LT 80-16-S	5 ml/h	4540004
Easypump® II LT 125-25-S	5 ml/h	4540006
Easypump® II LT 270-27-S	10 ml/h	4540008
Easypump® II LT 60-30-S	2 ml/h	4540010
Easypump® II LT 120-30-S	4 ml/h	4540012
Easypump® II LT 400-40-S	10 ml/h	4540014
Easypump® II LT 100-50-S	2 ml/h	4540016
Easypump® II LT 270-54-S	5 ml/h	4540018
Easypump® II LT 400-80-S	5 ml/h	4540022
Easypump® II LT 270-68-S	4 ml/h	4540026
Easypump® II LT 400-100-S	4 ml/h	4540028
Easypump® II LT 270-135-S	2 ml/h	4540032

## SHORT-TERM INFUSION

Designation	Color Code	Code No.* (REF)
Easypump® II ST 100-0.5-S	200 ml/h	4540040
Easypump® II ST 250-0.5-S	500 ml/h	4540042
Easypump® II ST 50-1-S	50 ml/h	4540044
Easypump® II ST 100-1-S	100 ml/h	4540046
Easypump® II ST 250-1-S	250 ml/h	4540048
Easypump® II ST 250-1.5-S	175 ml/h	4540050
Easypump® II ST 400-2-S	200 ml/h	4540052
Easypump® II ST 500-2-S	250 ml/h	4540054
Easypump® II ST 100-2-S	50 ml/h	4540056
Easypump® II ST 400-4-S	100 ml/h	4540058

All Easypump® II ST-LT versions are provided in boxes of 10 units.

\* suffix of article code may vary according to sales region

## SYSTEM AND CODING

Easypump® II offers a broad product portfolio of different versions. A clear coding system gives you guidance to identify the version that fits best your needs.

**LT** stands for "long-term" and means an infusion time between 12 hours and 135 hours.  
**ST** stands for "short-term" and means the range of pumps with an infusion time from 30 minutes to a maximum of 4 hours.

**30** stands for the nominal infusion time in hours. The following durations are the standard options:  
LT: 12-12.5-16-25-27-30-40-50-54-68-80-100-135 hours  
ST: 0.5-1-1.5-2-4 hours

Easypump® II LT 120-30-S

**120** indicates the nominal filling volume in ml. Easypump® II is available from 50 ml to 500 ml in the following volumes:  
LT: 60-80-100-120-125-270-400-500 ml  
ST: 50-100-250-400-500 ml

**S** means "soft shell", which describes a soft and flexible housing and body of the pump.



## COMPLEMENTARY PRODUCTS

Designation	Description	Units	Code No. (REF)
<b>Carry Pouches</b>			
Large	Comfortable transportation of large- and extra large-sized articles	1	4434455
Small	Comfortable transportation of small- and medium-sized articles	1	4434447
Insulated	Comfortable transportation of large- and extra large-sized articles which require a certain amount of temperature control	1	4438102