

Novosyn®

Absorbable braided suture
made of Polyglactin 910




Novosyn[®]

ABSORBABLE SUTURES

Description

Novosyn[®] is a mid-term absorbable braided and coated synthetic suture made of polyglactin 910. Novosyn[®] degrades by hydrolysis and provides a predictable and reliable absorption.




Novosyn[®]

Product overview

Polyglactin 910, braided, coated






Color	Violet or undyed
Origin	Synthetic
Sizes	USP 8/0 (metric 0.4) to USP 2 (metric 5)
Type of absorption	By hydrolysis
Mass absorption	From 56 to 70 days
Sterilisation	Ethylene oxide


	Thread length cm colour	USP (metric)							
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)
Pre-cut sutures									
5 threads of 70 cm	5 x 70 violet						B0058449	B0058450	
10 threads of 60 cm	10 x 60 violet					B0058567	B0058568		
6 threads of 45 cm	6 x 45 violet					B0058647	B0058649	B0058650	
12 threads of 45 cm	12 x 45 violet					B0058667	B0058669	B0058670	
1 thread of 150 cm	1 x 150 violet		B0058704	B0058705	B0058706	B0058707	B0058709	B0058708	
5 threads of 70 cm	5 x 70 undyed							B0059450	
6 threads of 45 cm	6 x 45 undyed					B0059647	B0059649		
12 threads of 45 cm	12 x 45 undyed					B0059667	B0059669	B0059670	
Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle									
	70 70 violet			B0068071		B0068072			
	70 violet		B0068073	B0068074	B0068075	B0068076	B0068077		
	150 violet					B0068119			

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Novosyn[®]

ABSORBABLE SUTURES



Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle with short cutting point									
 HRC 37 Easyslide	4 x 70 violet					B0068124			
1/2 circle round bodied needle									
 HR 37s Easyslide	150 violet						B0068125		
1/2 circle round bodied needle with short cutting point									
 HRC 37 Easyslide	90 violet						B0068143		
1/2 circle round bodied needle									
 HR 70 Easyslide	70 violet					B0068183			
1/2 circle round bodied needle with blunt point									
 HRN 50s Easyslide	90 violet						B0068229		
1/2 circle reverse cutting needle									
 HS 37s + HR 37s Easyslide	150 violet			B0068278					
 HS 48 Easyslide	150 violet					B0068303			




Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle with short cutting point									
 HRC 37ss 	90 violet				B0068344				
1/2 circle round bodied needle with trocar point									
 HRT 43  loop	150 violet				B0068346	B0068347			
1/2 circle round bodied needle									
 HR 40s  loop	150 violet					B0068350			
1/2 circle round bodied needle with short cutting point									
 HRC 37ss 	90 violet					B0068354			
1/2 circle round bodied needle									
 HR 40s  loop	150 violet				B0068355				
1/2 circle round bodied needle with trocar point									
 HRT 43s  loop	150 violet				B0068356	B0068357			

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Novosyn[®]

ABSORBABLE SUTURES






Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle with short cutting point									
 HRC 37ss	90 violet						B0068358		
Easyslide									
1/2 circle round bodied needle with trocar point									
 HRT 48	150 violet				B0068369	B0068370			
Easyslide									
1/2 circle round bodied needle									
 HR 60	150 violet					B0068380			
Easyslide									
1/2 circle round bodied needle with short cutting point									
 HRC 43ss	90 90 violet					B0068399	B0068404		
Easyslide									
3/8 circle round bodied needle									
 DR 48	90 violet					B0068437			
Easyslide									
1/2 circle reverse cutting needle									
 HS 48	90 violet						B0068463		
Easyslide									


Needle length mm	Thread length cm colour	USP (metric)								
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)	
1/2 circle round bodied needle with trocar point										
 HRT 43s 	90 violet						B0068473			
1/2 circle reverse cutting needle										
 HS 40 	90 violet						B0068488			
 HS 37 	90 violet						B0068498			
1/2 circle round bodied needle with trocar point										
 HRT 48 	90 150 violet					B0068520	B0068503			
 HRT 37s 	90 violet						B0068535			
1/2 circle round bodied needle										
 HR 40s 	120 90 violet					B0068555	B0068558			
 HR 48 	90 violet						B0068563			

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Novosyn[®]




ABSORBABLE SUTURES






Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle with short cutting point									
 HRC 48	90 violet						B0068575		
1/2 circle round bodied needle									
 HR 76	90 violet					B0068579			
1/2 circle reverse cutting needle									
 HS 40s	90 violet						B0068588		
1/2 circle round bodied needle									
 HR 37s	90 violet						B0068598		
1/2 circle reverse cutting needle									
 HS 37s	90 violet						B0068599		
1/2 circle round bodied needle with short cutting point									
 HRC 43s	150 violet				B0068601	B0068602			
3/8 circle round bodied needle with blunt point									
 DRN 64	90 violet				B0068607	B0068605			

Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle reverse cutting needle									
 HS 37s 	120 violet				B0068611				
3/8 circle round bodied needle with blunt point									
 DRN 64 	90 violet						B0068614		
1/2 circle reverse cutting needle									
 HS 37s 	120 violet					B0068632			
1/2 circle round bodied needle									
 HR 65 	90 violet					B0068639	B0068640		
1/2 circle round bodied needle with blunt point									
 HRN 45 	90 violet						B0068648		
1/2 circle round bodied needle with trocar point									
 HRT 37 	90 violet						B0068654		
1/2 circle round bodied needle									
 HR 37ss 	90 violet						B0068655		

Novosyn®




ABSORBABLE SUTURES





Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle									
	90 violet						B0068661		
1/2 circle round bodied needle with blunt point									
	90 violet					B0068674	B0068675		
1/2 circle round bodied needle									
	120 violet						B0068779		
	70 violet			B0068803	B0068804				
1/2 circle round bodied needle with short cutting point									
	150 violet					B0068837			
1/2 circle round bodied needle									
	90 violet						B0068848		

Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle with short cutting point									
	150 violet				B0068872	B0068873			
	90 violet						B0068897		
1/2 circle round bodied needle with trocar point									
	150 violet						B0068918		
1/2 circle round bodied needle									
	150 violet					B0068919			
	100 violet					B0068926			
	150 violet					B0068927	B0068928		

Novosyn®



ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)								
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)	
1/2 circle round bodied needle										
	240 violet					B0068929	B0068930			
	150 violet						B0068931			
	150 violet					B0068939	B0068940			
1/2 circle round bodied needle with short cutting point										
	90 violet				B0068950					
	90 violet						B0068983			
	150 violet				B0068986	B0068987	B0068988			
	90 violet						B0068992			

Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle									
	70 undyed			B0069064					
1/2 circle reverse cutting needle									
	90 undyed					B0069297			
	90 undyed					B0069364			
1/2 circle round bodied needle with trocar point									
	150 undyed					B0069370			
1/2 circle round bodied needle									
	140 undyed					B0069372			
	90 undyed					B0069558			
	90 undyed					B0069598			

Novosyn®

ABSORBABLE SUTURES


Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle with short cutting point									
 Easyslide	90 undyed				B0069886				
1/2 circle round bodied needle									
 Easyslide	90 violet						B0088002		
 Easyslide	4 x 70 90 violet						B0088003 B0088004		
1/2 circle round bodied needle with short cutting point									
 Easyslide	4 x 70 violet		B0088065						
1/2 circle round bodied needle									
 Easyslide	5 x 45 violet			B0088350					
 Easyslide	5 x 45 violet			B0088351					
 Easyslide	5 x 45 violet				B0088352				
 Easyslide	5 x 45 violet				B0088353	B0088354			

Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle									
 HR 40s take off	5 x 45 violet			B0088355					
 HR 22 take off	4 x 70 violet			B0088631					
 HR 37s take off	4 x 90 violet					B0088653			
 HR 65 take off	4 x 70 violet					B0088655	B0088656		
 HR 30 take off	4 x 70 violet			B0088661	B0088662				
 HR 40s take off	4 x 70 violet					B0088667	B0088668		
 HR 37 take off	4 x 70 violet			B0088670	B0088671				
 HR 65 take off	5 x 70 violet						B0088684		
 HR 48 take off	4 x 70 violet					B0088687	B0088688		

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Novosyn®

ABSORBABLE SUTURES


Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle									
 HR 26s Easyslide take off	4 x 70 violet				B0088691	B0088692	B0088693		
 HR 37s Easyslide take off	4 x 70 violet			B0088695		B0088697	B0088698		
1/2 circle reverse cutting needle									
 HS 37ss Easyslide take off	3 x 45 violet				B0088701				
 HS 37s Easyslide take off	4 x 45 violet					B0088705			
1/2 circle round bodied needle									
 HR 37ss Easyslide take off	4 x 70 violet					B0088707			
 HR 40s Easyslide take off	4 x 45 violet				B0088766	B0088767	B0088768		
 HR 43 Easyslide take off	4 x 45 violet					B0088781	B0088782		
 HR 48 Easyslide take off	4 x 45 violet			B0088785	B0088786	B0088787	B0088788		
 HR 26s Easyslide take off	4 x 45 violet				B0088791	B0088792	B0088793		

Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle round bodied needle									
 Easyslide take off	4 x 45 violet				B0088796	B0088797	B0088798		
 Easyslide take off	4 x 70 violet					B0088803			
 Easyslide take off	4 x 45 violet				B0088806	B0088807			
1/2 circle reverse cutting needle									
 Easyslide take off	4 x 45 violet				B0088906				
 Easyslide take off	3 x 45 4 x 45 violet					B0088911 B0088912			
1/2 circle round bodied needle									
 Easyslide take off	4 x 70 undyed				B0089696				
 Easyslide take off	4 x 45 undyed				B0089786				
 Easyslide take off	4 x 45 undyed				B0089796	B0089797			

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Novosyn®

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)							
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
1/2 circle reverse cutting needle									
	4 x 45 undyed					B0089912			
	Thread length cm colour	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)
Pre-cut sutures									
4 threads of 90 cm	4 x 90 violet				C0058064				
10 threads of 60 cm	10 x 60 violet	C0058180	C0058181	C0058182	C0058183				
1 thread of 140 cm	1 x 140 violet	C0058203	C0058204	C0058205	C0058207	C0058208	C0058209	C0058210	
2 threads of 70 cm	2 x 70 violet		C0058404	C0058405	C0058406	C0058407	C0058409	C0058410	
5 threads of 70 cm	5 x 70 violet		C0058444	C0058445	C0058446	C0058447			
5 threads of 90 cm	5 x 90 violet				C0058456	C0058457			
3 threads of 45 cm	3 x 45 violet	C0058623	C0058624	C0058625	C0058626	C0058627	C0058629	C0058630	
6 threads of 45 cm	6 x 45 violet	C0058643	C0058644	C0058645	C0058646				
12 threads of 45 cm	12 x 45 violet		C0058664	C0058665	C0058666				
1 thread of 140 cm	1 x 140 undyed			C0059205	C0059207	C0059208			
2 threads of 70 cm	2 x 70 undyed			C0059405			C0059409	C0059410	
5 threads of 70 cm	5 x 70 undyed			C0059445	C0059446	C0059447			
3 threads of 45 cm	3 x 45 undyed		C0059624	C0059625	C0059626	C0059627	C0059629	C0059630	
6 threads of 45 cm	6 x 45 undyed		C0059644	C0059645	C0059646				
12 threads of 45 cm	12 x 45 undyed		C0059664	C0059665	C0059666				

Needle length mm	Thread length cm colour	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
1/2 circle round bodied needle									
 HR 10 	70 70 violet		C0068002	C0068003					
 HR 22 	20 violet						C0068004		
 HR 13 	70 violet		C0068006	C0068007	C0068008	C0068009			
 HR 17 	70 violet		C0068011	C0068012	C0068013	C0068014	C0068015	C0068016	
 HR 26 	20 violet						C0068018		
 HR 26s 	45 violet						C0068019		
 HR 6 	45 violet	C0068020							
 HR 37s 	45 violet							C0068021	
3/8 circle reverse cutting needle with precision point									
 DSMP 24 	45 violet						C0068022		
1/2 circle round bodied needle with trocar point									
 HRT 40 	70 violet						C0068023	C0068024	
3/8 circle reverse cutting needle with precision point									
 DSMP 24 	45 violet					C0068025			

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Novosyn®

ABSORBABLE SUTURES


Needle		Thread	USP (metric)							
length mm		length cm colour	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
5/8 circle round bodied needle										
		70 violet					C0068027			
1/2 circle round bodied needle										
		70 violet			C0068028	C0068029				
		70 violet				C0068029N1				
		70 violet					C0068030			
		70 violet					C0068030N1			
		70 violet						C0068031		
		70 violet						C0068031N1		
		70 violet							C0068032	
		70 violet								
										
		70 violet								C0068033
										
Needle		Thread	USP (metric)							
length mm		length cm colour	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)
		70 violet							C0068034	
		70 violet		C0068035						
		70 violet				C0068036				
		70 violet				C0068036N1				
		70 violet					C0068037			
		70 violet					C0068037N1			
		70 violet						C0068038		
		70 violet							C0068038N1	
		70 violet								
										

Needle		Thread		USP (metric)								
length mm		length cm colour		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	
1/2 circle round bodied needle												
		70 violet	C0068039	C0068040								
		70 violet		C0068040N1								
		70 violet			C0068041							
		70 violet			C0068041N1							
		70 violet				C0068042						
		70 violet				C0068042N1						
		70 violet					C0068043					
		70 violet					C0068043N1					
		70 violet						C0068044				
		70 violet						C0068044N1				
			70 violet			C0068046				C0068045		
			70 violet			C0068046N1						
		70 violet				C0068047						
		70 violet				C0068047N1						
		70 violet					C0068048					
		70 violet					C0068048N1					
		70 violet						C0068049				
		70 violet						C0068049N1				

Novosyn®

ABSORBABLE SUTURES

Needle		Thread	USP (metric)								
length mm		length cm colour	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	
1/2 circle round bodied needle											
 HR 37		70 violet			C0068050						
		70 violet			C0068050N1						
		70 violet				C0068051					
		70 violet				C0068051N1					
		70 violet					C0068052				
		70 violet					C0068052N1				
		70 violet						C0068053			
		70 violet						C0068053N1			
		70 violet							C0068054		
											
	 HR 40s		70 violet				C0068055				
			70 violet				C0068055N1				
		70 violet					C0068056				
		70 violet					C0068056N1				
		70 violet						C0068057			
		70 violet						C0068057N1			
		70 violet			C0068059				C0068058		
											
 HR 48		70 violet				C0068060					
		70 violet				C0068060N1					
		70 violet					C0068061				

Needle length mm	Thread length cm colour	USP (metric)							
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)
1/2 circle round bodied needle									
 HR 48	70 violet					C0068061N1			
	70 violet						C0068062		
	70 violet						C0068062N1		
	70 violet							C0068063	
 HR 60	70 violet				C0068064	C0068065	C0068066		
 HR 26	120 violet				C0068067				
 HR 60	70 violet							C0068068	
 HR 22	20 violet			C0068069					
 HR 43s	70 violet							C0068070	
 HR 22	120 violet			C0068078					
 HR 26	20 violet			C0068079					
1/2 circle round bodied needle with short cutting point									
 HRC 17	70 violet						C0068081		

Novosyn®

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)							
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)
1/2 circle round bodied needle									
 HR 48	70 violet								
	70 violet					C0068061N1			
	70 violet						C0068062		
	70 violet						C0068062N1		
	70 violet							C0068063	
 HR 60	70 violet				C0068064	C0068065	C0068066		
 HR 26	120 violet				C0068067				
 HR 60	70 violet							C0068068	
 HR 22	20 violet			C0068069					
 HR 43s	70 violet							C0068070	
 HR 22	120 violet			C0068078					
 HR 26	20 violet			C0068079					
1/2 circle round bodied needle with short cutting point									
 HRC 17	70 violet						C0068081		

Needle length mm	Thread length cm colour	USP (metric)							
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)
1/2 circle reverse cutting needle									
 	HS 26ss 70 violet							C0068084	
 	HS 37ss 70 violet				C0068085	C0068086	C0068087		
1/2 circle round bodied needle									
 	HR 26 120 violet			C0068088					
 	 70 violet			C0068089					
	 70 violet			C0068089N1					
	 70 violet				C0068090				
	 70 violet				C0068090N1				
	 70 violet					C0068091			
	 70 violet					C0068091N1			
	 70 violet						C0068092		
	 70 violet						C0068092N1		
	 70 violet							C0068093	

Novosyn®

ABSORBABLE SUTURES






Needle length mm	Thread length cm colour	USP (metric)							
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)
1/2 circle round bodied needle									
 HR 37s	70 violet			C0068094					
	 70 violet			C0068094N1					
	70 violet				C0068095				
	 70 violet				C0068095N1				
	70 violet					C0068096			
	 70 violet					C0068096N1			
	70 violet						C0068097		
	 70 violet						C0068097N1		
	 70 violet							C0068098	
 HR 48	70 violet			C0068099					
									
 HR 26ss	70 violet						C0068100		
									
Needle length mm	Thread length cm colour	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
3/8 circle round bodied needle									
 DR 10	45 violet	C0068101	C0068102	C0068103					
									
 DR 12	60 violet				C0068105				
									
 DR 15	70 violet				C0068107				
									
 DR 22 b	70 violet				C0068108				
									

Needle length mm	Thread length cm colour	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
1/2 circle round bodied needle									
 	45 violet						C0068109		
 	70 violet						C0068110	C0068111	
 	70 violet				C0068113	C0068114	C0068115	C0068116	
3/8 circle round bodied needle									
 	70 violet					C0068117			
1/2 circle round bodied needle									
 	70 violet								C0068118
 	70 violet					C0068120			
 	70 violet						C0068121	C0068122	C0068123
Needle length mm	Thread length cm colour	USP (metric)							
		8/0 (0.4)	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)
3/8 circle round bodied needle									
 	30 violet	C0068133							
 	45 violet				C0068134				
 	45 violet					C0068137			
 	45 violet					C0068138			

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Novosyn[®]

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)							
		8/0 (0.4)	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)
1/2 circle round bodied needle with short cutting point									
	90 violet						C0068139	C0068140	
	90 violet							C0068140N1	
	90 violet								C0068141
	90 violet								C0068141N1
Needle length mm	Thread length cm colour	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
	90 violet								C0068142
	90 violet								C0068142N1
1/2 circle round bodied needle									
	70 violet			C0068144	C0068147				
	30 violet	C0068149							
3/8 circle round bodied needle									
	70 violet				C0068150	C0068151			
1/2 circle round bodied needle									
	70 violet					C0068152			
Ski Needle round bodied needle									
	45 violet					C0068154			
3/8 circle round bodied needle									
	70 violet			C0068155					


Needle length mm	Thread length cm colour	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
3/8 circle round bodied needle									
 DR 22	70 violet				C0068159	C0068160			
	70 violet					C0068160N1			
	70 violet						C0068161		
									
3/8 circle reverse cutting needle with precision point									
 2xDSMP 13	45 violet			C0068162					
									
1/2 circle round bodied needle									
 2xHR 17b	70 violet			C0068163					
									
3/8 circle round bodied needle									
 DR 26	70 violet					C0068164	C0068165		
	70 violet						C0068165N1		
									
Needle length mm	Thread length cm colour	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
 DR 32	70 violet								C0068168
									
1/2 circle round bodied needle									
 HR 43s	70 violet							C0068169	
									
3/8 circle round bodied needle									
 DR 32	70 violet				C0068170	C0068171			
	70 violet					C0068171N1			
	70 violet						C0068172	C0068173	
									

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Novosyn®

ABSORBABLE SUTURES



Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle									
 HR 40	70 violet						C0068175		
	70 violet						C0068175N1		
 2xHR 17b	70 violet			C0068176					
3/8 circle round bodied needle									
 DR 26	120 violet				C0068179				
1/2 circle round bodied needle									
 HR 22b	70 violet			C0068181					
Straight round bodied needle									
 GR 51	90 violet					C0068185			
1/2 circle round bodied needle									
 HR 22b	70 violet				C0068186				
	70 violet				C0068186N1				
1/2 circle round bodied needle with short cutting point									
 HRC 37	150 violet					C0068187			
1/2 circle round bodied needle									
 HR 22b	70 violet					C0068189			
 2xHR 22b	70			C0068190					
	90			C0068191					
	70 violet				C0068192				

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
Hook round bodied needle with short cutting point									
	120 violet					C0068193	C0068194		
1/2 circle round bodied needle									
	70 violet			C0068196	C0068198				
	70 violet				C0068198N1				
3/8 circle reverse cutting needle									
	45 violet	C0068202							
1/2 circle reverse cutting needle									
	45 violet		C0068203	C0068204					
	70 violet							C0068205	
3/8 circle reverse cutting needle									
	45 violet	C0068206	C0068207	C0068208					
1/2 circle round bodied needle									
	70 violet					C0068209			
	70 violet					C0068209N1			
	70 violet				C0068210				
Straight reverse cutting needle									
	45 violet			C0068211					

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Novosyn®

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
3/8 circle reverse cutting needle									
 DS 16	45 violet		C0068212	C0068213	C0068214				
1/2 circle round bodied needle									
 2xHR 26b	90 70 violet				C0068216	C0068217			
3/8 circle reverse cutting needle									
 DS 19	45 violet		C0068219	C0068220					
	45 violet			C0068220N1					
	45 violet				C0068221				
	45 violet				C0068221N1				
	90 violet			C0068224					
1/2 circle reverse cutting needle									
 HS 21	45 violet			C0068225	C0068226				
 HS 55	70 violet								C0068227
1/2 circle round bodied needle with blunt point									
 HRN 50s	90 violet							C0068228	
1/2 circle reverse cutting needle									
 HS 23	45 violet				C0068230				

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
3/8 circle reverse cutting needle									
 DS 24 	70 violet		C0068232						
1/2 circle reverse cutting needle									
 HS 15 	70 violet		C0068233						
3/8 circle reverse cutting needle									
 DS 24 	70 violet			C0068234					
	 70 violet			C0068234N1					
	70 violet				C0068235				
	 70 violet				C0068235N1				
	70 violet					C0068236			
	 70 violet					C0068236N1			
	70 violet						C0068237		
	 70 violet						C0068237N1		
	90 violet						C0068239	C0068240	
1/2 circle reverse cutting needle									
 HS 26 	70 violet				C0068241				
	 70 violet				C0068241N1				
	70 violet					C0068242			
	 70 violet					C0068242N1			
	70 violet						C0068243	C0068244	

Novosyn®

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
3/8 circle reverse cutting needle									
 DS 30	70 violet				C0068246				C0068245
	 70 violet				C0068246N1				
	70 violet					C0068247			
	 70 violet					C0068247N1			
	70 violet						C0068248		
	 70 violet						C0068248N1		
	 70 violet							C0068249	
1/2 circle reverse cutting needle									
 HS 15	70 violet			C0068250					
 HS 30	70 violet					C0068251			
	 70 violet					C0068251N1			
	70 violet						C0068252		
	 70 violet						C0068252N1		
	 70 violet							C0068253	
3/8 circle reverse cutting needle									
 DS 35	70 violet				C0068255	C0068256	C0068257	C0068258	C0068259
1/2 circle reverse cutting needle									
 HS 48	70 violet						C0068261	C0068262	C0068263
Straight reverse cutting needle									
 GS 60	70 violet			C0068264	C0068265	C0068266			

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle reverse cutting needle									
 	70 violet					C0068267	C0068268		
	70 violet						C0068268N1		
	70 violet							C0068269	
	70 violet							C0068269N1	
Straight reverse cutting needle									
 	70 violet				C0068272	C0068274	C0068275		
Straight round bodied needle									
 	20 violet			C0068281					
3/8 circle reverse cutting needle									
 	70 violet						C0068282	C0068283	
Straight round bodied needle									
 	20 70 violet				C0068284 C0068285				
Ski Needle round bodied needle									
 	20 violet					C0068287			
Straight round bodied needle									
 	20 violet			C0068288					

Novosyn®

ABSORBABLE SUTURES

Needle		Thread	USP (metric)							
length mm		length cm colour	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle reverse cutting needle										
 HS 26s	 	70 violet					C0068290			
		70 violet					C0068290N1			
		70 violet						C0068291		
		70 violet						C0068291N1		
		70 violet							C0068292	
		70 violet								
 HS 37s	  	70 violet				C0068293	C0068294			
		70 violet					C0068294N1			
		70 violet						C0068296		
		70 violet						C0068296N1		
		70 violet							C0068297	
		70 violet							C0068297N1	
		70 violet								C0068298
		70 violet								
3/8 circle round bodied needle with short cutting point										
 DRC 12		70 violet	C0068300	C0068301	C0068302					
3/8 circle round bodied needle with trocar point										
 DRT 12		70 violet	C0068306	C0068307						
 DRT 32		70 violet							C0068309	
1/2 circle round bodied needle with short cutting point										
 HRC 22		120 violet				C0068310				
 HRC 17		70 violet		C0068312	C0068313	C0068314	C0068315	C0068316		






Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
3/8 circle round bodied needle with short cutting point									
  DRC 26	70 violet					C0068317			
3/8 circle round bodied needle with trocar point									
  DRT 18	70 violet		C0068319	C0068320	C0068321				
3/8 circle round bodied needle with short cutting point									
  DRC 32	70 70 violet					C0068322	C0068324		
1/2 circle round bodied needle with short cutting point									
  HRC 22	70 violet			C0068325	C0068326				
  HRC 37ss	70 violet					C0068329			
  HRC 22	70 violet					C0068331	C0068332	C0068333	
  HRC 26	20 70 violet				C0068334 C0068335				
	70 violet				C0068335N1				
	70 violet					C0068336			
	70 violet					C0068336N1			
	20 violet					C0068337			

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Novosyn®

ABSORBABLE SUTURES


Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle with short cutting point									
	70 70 90 violet					C0068342	C0068338	C0068340	
									
	70 violet						C0068343		
									
	90 violet						C0068344N1		
									
	120 violet					C0068345			
									
1/2 circle round bodied needle with trocar point									
	70 violet					C0068349			
									
	70 violet					C0068352			C0068353
									
1/2 circle round bodied needle with short cutting point									
	90 violet							C0068354N1	
									




Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle with short cutting point										
	70 violet							C0068359		
	70 violet				C0068360	C0068361				
	70 violet					C0068361N1				
	70 violet						C0068362	C0068363		
1/2 circle round bodied needle										
	120 violet				C0068365					
1/2 circle round bodied needle with short cutting point										
	70 violet								C0068368	
	70 violet					C0068372	C0068371			
	70 violet					C0068372N1				
	70 violet				C0068373			C0068374		
	90 violet							C0068375		
3/8 circle round bodied needle with trocar point										
	70 violet			C0068376	C0068377					
3/8 circle reverse cutting needle										
	70 violet					C0068379				

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Novosyn®

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle									
 HR 40 Easyslide	70 violet					C0068381			
3/8 circle round bodied needle with trocar point									
 DRT 26 Easyslide	70 violet				C0068382	C0068383			
1/2 circle round bodied needle with trocar point									
 HRT 37s Easyslide	90 violet				C0068384				
3/8 circle round bodied needle with trocar point									
 DRT 32 Easyslide	70 violet					C0068386	C0068387		
1/2 circle round bodied needle with trocar point									
 HRT 37 Easyslide	70 violet					C0068392	C0068393	C0068391	
 HRT 37s Easyslide	70 violet					C0068395	C0068396	C0068397	
5/8 circle reverse cutting needle									
 FS 25 Easyslide	70 violet				C0068398				
1/2 circle round bodied needle with trocar point									
 HRT 48 Easyslide	70 violet					C0068400	C0068401	C0068402	C0068403

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
Straight round bodied needle									
 GR 51	70 violet				C0068405				
1/2 circle round bodied needle with short cutting point									
 HRC 43ss	70 violet							C0068406	
Needle length mm	Thread length cm colour	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
3/8 circle reverse cutting needle with precision point									
 DSMP 7	70 violet	C0068407							
Hook round bodied needle									
 JR 36ss	70 violet							C0068408	
3/8 circle round bodied needle									
 DR 40	70 violet								C0068409
Hook round bodied needle									
 JR 36ss	70 violet						C0068410		
1/2 circle round bodied needle with trocar point									
 HRT 40s	70 violet							C0068413	

Novosyn®


ABSORBABLE SUTURES


Needle		Thread	USP (metric)							
length mm		length cm colour	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
5/8 circle round bodied needle										
 	FR 26	70 violet					C0068414	C0068415	C0068416	C0068417
3/8 circle reverse cutting needle										
 	DS 19	70 violet			C0068419	C0068420				
		 70 violet				C0068420N1				
		70 violet					C0068421			
		 70 violet					C0068421N1			
		70 violet						C0068422		
		 70 violet						C0068422N1		
5/8 circle round bodied needle										
 	FR 36	70 violet					C0068424	C0068425	C0068426	C0068427
 	FR 26s	70 violet								C0068428
Needle		Thread	USP (metric)							
length mm		length cm colour	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle reverse cutting needle										
 	HS 23	70 violet				C0068430				C0068429
		 70 violet				C0068430N1				
		70 violet					C0068431			
		 70 violet					C0068431N1			
		70 violet						C0068432		
		 70 violet								
3/8 circle round bodied needle										
 	DR 8	30 violet	C0068433							

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
3/8 circle round bodied needle									
 DR 12 	70 violet	C0068434	C0068435	C0068436					
1/2 circle round bodied needle									
 HR 22 	90 violet			C0068438					
1/2 circle reverse cutting needle									
 HS 18 	70 violet		C0068439						
 HS 26s 	90 violet					C0068440			
1/2 circle round bodied needle with short cutting point									
 HRC 26s 	70 violet					C0068441			
1/2 circle reverse cutting needle									
 HS 26 	90 violet					C0068442			
1/2 circle round bodied needle									
 HR 22 	70 violet								C0068443
1/2 circle round bodied needle with short cutting point									
 HRC 26s 	70 violet						C0068444		

Novosyn®




ABSORBABLE SUTURES


Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle									
 	90 violet				C0068446	C0068447			
	90 violet					C0068447N1			
	90 violet						C0068448		
	90 violet						C0068448N1		
 	45	C0068449	C0068450	C0068451					
	90 violet			C0068452					
3/8 circle reverse cutting needle									
 	70 violet							C0068453	
1/2 circle round bodied needle									
 	70 violet							C0068457	
	70 violet							C0068457N1	
1/2 circle reverse cutting needle									
 	90 violet				C0068458				
3/8 circle round bodied needle with short cutting point									
 	70 violet			C0068460					
1/2 circle reverse cutting needle									
 	90 violet						C0068461	C0068462	
	90 violet							C0068462N1	

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle reverse cutting needle									
 	HS 48f 90 violet							C0068464	
Straight reverse cutting needle									
 	GS 60 45 violet				C0068465				
1/2 circle round bodied needle									
  	HR 40s 90 violet				C0068466	C0068467			
	90 violet					C0068467N1			
3/8 circle round bodied needle with short cutting point									
 	DRC 18 70 violet				C0068468				
1/2 circle round bodied needle with trocar point									
 	HRT 43s 90 violet							C0068470	
 	HRT 43 90 violet						C0068471	C0068472	
 	HRT 43s 70 violet							C0068474	
Hook round bodied needle									
 	JR 30 70 violet					C0068475			

Novosyn®

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
Hook round bodied needle with short cutting point									
 	70 violet						C0068477	C0068476	
Hook round bodied needle									
 	70 violet							C0068478	
Hook round bodied needle with short cutting point									
 	45 violet						C0068479		
1/2 circle reverse cutting needle									
 	90 violet						C0068480		
Hook round bodied needle									
 	70 violet							C0068481	
Hook round bodied needle with short cutting point									
 	60 violet					C0068483			
1/2 circle round bodied needle with trocar point									
 	70 violet							C0068484	

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle reverse cutting needle										
  	90 violet					C0068485	C0068486			
	90 violet						C0068486N1			
	90 violet							C0068487		
	90 violet							C0068487N1		
	90 violet									
5/8 circle round bodied needle										
 	70 violet					C0068489	C0068490	C0068491		
1/2 circle round bodied needle with trocar point										
 	70 violet								C0068492	
1/2 circle reverse cutting needle										
 	90 violet							C0068493		
    	90 violet				C0068494	C0068495				
	90 violet					C0068495N1				
	90 violet						C0068496			
	90 violet						C0068496N1			
	90 violet							C0068497		
	90 violet							C0068497N1		
	90 violet									
1/2 circle round bodied needle with trocar point										
 	90 violet				C0068504	C0068500	C0068501	C0068502		

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Novosyn®

ABSORBABLE SUTURES



Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
5/8 circle round bodied needle with short cutting point										
 FRC 22 	90 violet				C0068505					
3/8 circle reverse cutting needle										
 DS 12 	70 violet	C0068506	C0068507	C0068508						
 DS 16 	90 70 violet		C0068511 C0068512	C0068513	C0068514					
1/2 circle round bodied needle with protection point										
 HRB 40s 	90 violet							C0068516		
1/2 circle round bodied needle with blunt point										
 HRN 35 	70 violet					C0068517	C0068518			
1/2 circle reverse cutting needle										
 HS 18 	70 violet			C0068522						
5/8 circle round bodied needle with short cutting point										
 FRC 22 	90 violet			C0068523						
1/2 circle reverse cutting needle										
 HS 23 	90 violet					C0068524				
1/2 circle round bodied needle with trocar point										
 HRT 40 	70 violet							C0068525		

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle									
 	90 violet				C0068526				
3/8 circle reverse cutting needle with precision point									
 	70 violet		C0068528						
 	45 violet	C0068530		C0068531	C0068529				
1/2 circle round bodied needle with trocar point									
 	90 violet					C0068532	C0068533	C0068534	
3/8 circle reverse cutting needle									
 	45 violet			C0068536					
3/8 circle reverse cutting needle with precision point									
 	45 violet		C0068538						
1/2 circle round bodied needle									
	90 violet			C0068540	C0068541				
	 90 violet				C0068541N1				
	90 violet					C0068542			
	 90 violet					C0068542N1			
	90 violet						C0068543		
	 90 violet						C0068543N1		
	 90 violet							C0068544	C0068545

Novosyn®

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle									
 	70 violet					C0068547	C0068548		
	70 violet						C0068548N1		
	70 violet							C0068549	
	70 violet							C0068549N1	
 	90 violet				C0068550				
	90 violet				C0068550N1				
	90 violet					C0068551			
	90 violet					C0068551N1			
	90 violet						C0068552		
	90 violet						C0068552N1		
	90 violet							C0068553	
	90 violet							C0068553N1	
 	90 violet						C0068556		
	90 violet						C0068556N1		
	90 violet							C0068557	
	90 violet							C0068557N1	
 	70 violet							C0068559	

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle										
 HR 48	90 violet					C0068560				
	90 violet					C0068560N1				
	90 violet						C0068561			
	90 violet						C0068561N1			
	90 violet							C0068562		
	90 violet							C0068562N1		
3/8 circle round bodied needle										
 DR 32	90 violet					C0068566	C0068567	C0068568		
1/2 circle round bodied needle										
 HR 40	90 violet						C0068569			
	90 violet						C0068569N1			
1/2 circle round bodied needle with short cutting point										
 HRC 48	90 violet					C0068570				
1/2 circle round bodied needle										
 HR 40	90 violet							C0068571		
	90 violet							C0068571N1		
1/2 circle round bodied needle with short cutting point										
 HRC 48	90 violet						C0068572	C0068573	C0068574	

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Novosyn®

ABSORBABLE SUTURES


Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle										
  	90 violet						C0068576	C0068577		
	90 violet							C0068577N1		
 	90 violet							C0068580		
1/2 circle reverse cutting needle										
  	90 violet						C0068582			
	90 violet						C0068582N1			
	90 violet							C0068583		
	90 violet							C0068583N1		
	90 violet									
1/2 circle round bodied needle										
 	150 violet					C0068584				
1/2 circle reverse cutting needle										
  	90 violet					C0068585	C0068586			
	90 violet						C0068586N1			
	90 violet							C0068587		
	90 violet							C0068587N1		
	90 violet									

Needle	Thread	USP (metric)							
length mm	length cm colour	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle									
	90 violet				C0068589	C0068590			
	90 violet					C0068590N1			
	90 violet						C0068591		
	90 violet						C0068591N1		
	90 violet							C0068592	
	90 violet							C0068592N1	
	90 violet								
1/2 circle reverse cutting needle									
	90 violet				C0068594				
1/2 circle round bodied needle									
	90 violet					C0068595			
	90 violet					C0068595N1			
1/2 circle reverse cutting needle									
	90 violet					C0068596			
	90 violet					C0068596N1			
1/2 circle round bodied needle									
	90 violet							C0068597	
	90 violet							C0068597N1	
1/2 circle reverse cutting needle									
	120 violet				C0068600				

Novosyn®


ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle reverse cutting needle										
 	120 violet					C0068603				
1/2 circle round bodied needle with blunt point										
 	70 violet					C0068604				
 	70 violet							C0068608		
1/2 circle reverse cutting needle										
 	120 violet						C0068610			
1/2 circle round bodied needle										
 	90 violet								C0068613	
1/2 circle round bodied needle with blunt point										
 	90 violet					C0068617	C0068618	C0068619		
3/8 circle reverse cutting needle with precision point										
 	45 violet	C0068622	C0068623	C0068624						
1/2 circle round bodied needle with short cutting point										
 	70 violet					C0068625				

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
3/8 circle reverse cutting needle with precision point									
	45 violet		C0068626						
1/2 circle round bodied needle with short cutting point									
	70 violet							C0068627	C0068628
Ski Needle round bodied needle									
	30 violet					C0068629			
3/8 circle reverse cutting needle with precision point									
	45 violet			C0068630	C0068631				
1/2 circle round bodied needle									
	90 violet							C0068633	
1/2 circle round bodied needle with blunt point									
	70 90 violet						C0068634 C0068636	C0068635 C0068637	
	70 violet				C0068638				
1/2 circle reverse cutting needle									
	70 violet				C0068641				

Novosyn®

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle										
	90 violet									C0068642
	90 violet								C0068645	
	90 violet								C0068646	
1/2 circle round bodied needle with blunt point										
	90 violet								C0068647	
1/2 circle round bodied needle with trocar point										
	90 violet				C0068650	C0068651	C0068652	C0068653		
3/8 circle reverse cutting needle with precision point										
	70 violet				C0068659	C0068660				
1/2 circle round bodied needle with blunt point										
	90 violet						C0068663	C0068664		
3/8 circle reverse cutting needle with precision point										
	70 violet				C0068669	C0068670	C0068671			

Needle length mm	Thread length cm colour	USP (metric)							
		8/0 (0.4)	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)
1/2 circle reverse cutting needle with precision point									
 2xHSMP 9 	35 violet	C0068673							
3/8 circle reverse cutting needle with precision point									
 DSMP 11 	45 70 violet		C0068684	C0068685 C0068688					
 DSMP 35 	70 violet							C0068690	
1/2 circle reverse cutting needle									
 HS 48 	70 violet							C0068691	
1/2 circle round bodied needle with short cutting point									
 HRC 26 	70 violet					C0068694			
Needle length mm	Thread length cm colour	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
Straight round bodied needle									
 GR 19 	70 violet								C0068695
3/8 circle reverse cutting needle									
 2xDS 12 	45 violet		C0068706						
3/8 circle reverse cutting needle with precision point									
 2xDSMP 7 	45 violet		C0068708						

Novosyn[®]

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
1/2 circle round bodied needle with short cutting point									
	70 violet							C0068716	
	70 violet							C0068721	
3/8 circle reverse cutting needle with precision point									
	70 violet			C0068723					
	45 violet		C0068726						
1/4 circle lancet needle with micro-point									
	45 violet	C0068732							
Ski Needle round bodied needle									
	120 violet					C0068744			
	40 violet					C0068745			
	20 violet					C0068746			
1/2 circle round bodied needle									
	70 violet						C0068747		
Ski Needle round bodied needle									
	20 120 violet						C0068749 C0068751		
	20 120 violet				C0068752 C0068755	C0068753 C0068756	C0068754 C0068757		

Needle length mm	Thread length cm colour	USP (metric)								
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	
1/2 circle reverse cutting needle with precision point										
2xHSMP 9 	45 violet		C0068758							
	1/2 circle round bodied needle									
	HR 20 	90 70 violet				C0068765	C0068763			
		70 violet				C0068765N1				
		70 violet					C0068766			
		70 violet					C0068766N1			
70 violet							C0068767			
3/8 circle reverse cutting needle										
2xDS 12 	45 violet			C0068769						
	Needle									
length mm	Thread length cm colour	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle reverse cutting needle										
HS 40s 	70 violet								C0068770	
	5/8 circle round bodied needle									
FR 26 	90 violet					C0068775				
	Ski Needle round bodied needle									
SKR 26 	120 violet						C0068781			
	3/8 circle reverse cutting needle with precision point									
DSMP 16 	70 violet		C0068782	C0068783	C0068784					

Novosyn®


ABSORBABLE SUTURES


Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle with short cutting point										
 HRC 20s 	70 violet					C0068785				
1/2 circle round bodied needle										
 HR 26ss 	70 violet						C0068791		C0068793	
5/8 circle round bodied needle with short cutting point										
 FRC 10 	30 violet	C0068798								
1/2 circle reverse cutting needle										
 HS 40 	70 violet						C0068805	C0068806		
Progressive curved round bodied needle with short cutting point										
 PRC 24 	120 violet					C0068807				
5/8 circle round bodied needle										
 FR 26 	90 violet						C0068810			
 2xFR 26 	70 90 violet			C0068813	C0068818	C0068815 C0068820	C0068817			
3/8 circle reverse cutting needle with precision point										
 DSMP 19 	70 violet			C0068823	C0068824					
	70 violet				C0068824N1					
3/8 circle round bodied needle with short cutting point										
 DRC 22 	70 violet				C0068827					

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
5/8 circle round bodied needle									
	90 violet						C0068828		
									
3/8 circle round bodied needle with short cutting point									
	20 violet					C0068830			
									
	20 violet				C0068833				
									
1/2 circle reverse cutting needle									
	70 violet								C0068840
									
1/2 circle round bodied needle with short cutting point									
	70 violet					C0068842	C0068843	C0068844	
									
1/2 circle round bodied needle									
	90 violet							C0068847	
	90 violet							C0068847N1	
									
1/2 circle round bodied needle with short cutting point									
	140 violet				C0068852				
									
1/2 circle round bodied needle									
	90 violet			C0068853	C0068854	C0068855			
									

Novosyn[®]

ABSORBABLE SUTURES




Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
Hook round bodied needle with short cutting point										
	90 violet							C0068857		
5/8 circle round bodied needle										
	90 violet							C0068860		
	90 violet							C0068861		
3/8 circle reverse cutting needle with precision point										
	70 violet				C0068864	C0068865				
	70 violet					C0068865N1				
5/8 circle round bodied needle										
	90 violet						C0068866	C0068867		
1/2 circle round bodied needle										
	90 violet		C0068869	C0068870						
	90 violet				C0068874	C0068875				
	90 violet				C0068877					
	90 violet				C0068877N1					
	120 90 violet					C0068880 C0068881				

Needle length mm	Thread length cm colour	6/0 (0.7)	5/0 (1)	4/0 (1.5)	USP (metric)		0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle									
 HR 17f 	70 violet			C0068888					
1/2 circle round bodied needle with short cutting point									
 HRC 37s 	90 violet				C0068893				
	 90 violet				C0068893N1				
	90 violet					C0068894			
	 90 violet					C0068894N1			
	90 violet						C0068895		
	 90 violet						C0068895N1		
	90 violet							C0068896	
	 90 violet							C0068896N1	
1/2 circle round bodied needle									
 2xHR 26 	90 violet				C0068901	C0068902			
 2xHR 17 	70 violet				C0068904				
 2xHR 26 	70 violet						C0068907		
1/2 circle round bodied needle with short cutting point									
 HRC 40 	70 violet					C0068908	C0068909		

Novosyn®





ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
Straight reverse cutting needle									
 2xGS 13 	45 violet			C0068920					
1/2 circle round bodied needle									
 HR 22s 	70 violet			C0068932	C0068934				
 2xHR 22s 	70 violet						C0068937		
Straight round bodied needle									
 2xGR 19 	35 violet			C0068943					
 GR 19 	35 violet			C0068944					
1/2 circle round bodied needle									
 HR 43 	70 violet						C0068946	C0068947	
	70 violet							C0068947N1	
	70 violet								C0068948
1/2 circle round bodied needle with short cutting point									
 HRC 40 	90 violet					C0068949		C0068951	
1/2 circle round bodied needle									
 HR 37s 	90 violet			C0068953	C0068955		C0068956		
	90 violet						C0068956N1		
1/2 circle reverse cutting needle									
 HS 26 	120 violet					C0068957			

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle with short cutting point									
 HRC 40s	70 violet					C0068958	C0068959	C0068960	
	70 violet							C0068960N1	
	70 violet								C0068961
 HRC 48	70 violet					C0068962	C0068963	C0068964	C0068965
 HRC 43	70 violet					C0068971	C0068972		
1/2 circle reverse cutting needle									
 HS 26ss	70 violet						C0068975		
1/2 circle round bodied needle with short cutting point									
 HRC 43	90 violet					C0068976	C0068977	C0068978	
 HRC 43s	90 violet							C0068982	
Straight round bodied needle									
 GR 65	90 violet					C0068994			

Novosyn[®]

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle with trocar point									
	90 violet					C0068997			
5/8 circle round bodied needle with short cutting point									
	70 undyed			C0069003					
1/2 circle round bodied needle with short cutting point									
	90 undyed					C0069004			
1/2 circle round bodied needle									
	70 undyed		C0069007	C0069008					
1/2 circle round bodied needle with short cutting point									
	70 undyed					C0069010			
1/2 circle round bodied needle									
	70 undyed	C0069011	C0069012	C0069013	C0069014	C0069015			
	70 undyed			C0069029	C0069030	C0069031		C0069033	

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle									
 HR 26 	70 undyed			C0069040	C0069041				
	70 undyed				C0069041N1				
	70 undyed					C0069042			
	70 undyed					C0069042N1			
	70 undyed						C0069043		
 HR 30 	70 undyed					C0069047	C0069048		
 HR 37 	70 undyed				C0069050	C0069051	C0069052		
 HR 40s 	70 undyed					C0069055		C0069057	
 HR 48 	70 undyed				C0069059				
1/2 circle reverse cutting needle									
 HS 26ss 	70 undyed							C0069083	
 HS 26s 	70 undyed						C0069084		
 HS 37ss 	70 undyed					C0069085	C0069088	C0069087	

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Novosyn®

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle									
 HR 26s 	70 undyed				C0069089	C0069090	C0069091	C0069092	
1/2 circle reverse cutting needle									
 HS 40s 	70 undyed							C0069093	
1/2 circle round bodied needle									
 HR 37s 	70 undyed				C0069099	C0069095	C0069096	C0069097	C0069098
Progressive curved round bodied needle									
 PR 16 	45 undyed			C0069135					
3/8 circle reverse cutting needle									
 DS 30 	90 undyed						C0069148		
1/2 circle reverse cutting needle									
 HS 37s 	90 70 undyed					C0069187 C0069195	C0069196	C0069197	
 HS 15 	45 undyed		C0069203	C0069204					
3/8 circle reverse cutting needle									
 DS 12 	45 undyed	C0069206	C0069207	C0069208					
 DS 16 	45 undyed		C0069212	C0069213	C0069214				

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
3/8 circle reverse cutting needle									
 DS 19f 	45 undyed			C0069217					
	 DS 19 	45 undyed		C0069219	C0069220				
45 undyed				C0069220N1					
45 undyed					C0069221				
45 undyed					C0069221N1				
90 undyed				C0069223					
1/2 circle reverse cutting needle									
 HS 23 	45 undyed				C0069230				
3/8 circle reverse cutting needle									
 DS 24 	90 70 undyed			C0069234	C0069235	C0069232			
	70 undyed				C0069235N1				
	70 undyed					C0069236			
	70 undyed					C0069236N1			
	70 undyed						C0069237		
	 DS 30 	70 undyed				C0069246	C0069247		C0069244
70 undyed						C0069247N1			
70 undyed							C0069248		
 DS 35 	70 undyed					C0069256			
Straight reverse cutting needle									
 GS 60 	70 undyed				C0069265				

C Box of 36 units

B Box of 24 units



G Box of 12 units

M Box of 6 units

Novosyn®

ABSORBABLE SUTURES


Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
Straight reverse cutting needle										
 GS 60f 	70 undyed				C0069266					
 GS 60 	70 undyed					C0069267				
 GS 51 	70 undyed			C0069271	C0069272	C0069273				
 GS 60 	70 undyed			C0069280						
3/8 circle reverse cutting needle										
 DS 48 	70 undyed								C0069282	
1/2 circle reverse cutting needle										
 HS 26s 	70 undyed					C0069290		C0069293		
 HS 37 	90 undyed					C0069294	C0069295	C0069296		
5/8 circle round bodied needle with short cutting point										
 FRC 22 	90 undyed			C0069299	C0069298					
3/8 circle reverse cutting needle										
 DS 16 	70 undyed				C0069302					
1/2 circle round bodied needle with short cutting point										
 HRC 17 	45 undyed		C0069309							
3/8 circle round bodied needle with trocar point										
 DRT 18 	70 undyed			C0069320	C0069321					

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle reverse cutting needle									
 HS 18 	70 undyed			C0069322					
1/2 circle round bodied needle with short cutting point									
 HRC 22 	70 undyed				C0069324				
3/8 circle reverse cutting needle with precision point									
 DSMP 16 	70 undyed		C0069325	C0069326	C0069327				
3/8 circle reverse cutting needle									
 DS 12 	70 undyed	C0069330	C0069331						
1/2 circle round bodied needle with short cutting point									
 HRC 22 	70 undyed			C0069335					
 HRC 17 	70 undyed		C0069339	C0069340	C0069341				
1/2 circle reverse cutting needle									
 HS 48 	90 undyed						C0069362	C0069363	
1/2 circle round bodied needle with trocar point									
 HRT 43s 	90 undyed					C0069378			
 HRT 37 	70 undyed					C0069392			

Novosyn®

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle with trocar point										
 HRT 40 ⊕ Easyslide	70 undyed					C0069400				
1/2 circle round bodied needle										
 HR 13 ● Easyslide	45 undyed	C0069405	C0069407	C0069408						
5/8 circle round bodied needle										
 FR 26 ● Easyslide	70 undyed					C0069415				
3/8 circle reverse cutting needle										
 DS 19 ▼ Easyslide	70 undyed		C0069419	C0069420						
	70 undyed			C0069420N1						
	70 undyed				C0069421					
	70 undyed				C0069421N1					
	70 undyed					C0069422				
	70 undyed									
1/2 circle reverse cutting needle										
 HS 21 ▼ Easyslide	70 undyed				C0069425					
 HS 23 ▼ Easyslide	70 undyed				C0069430	C0069431	C0069432	C0069433		
 HS 23ss ▼ Easyslide	70 undyed							C0069438		
 HS 26 ▼ Easyslide	70 undyed				C0069440	C0069441				

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle with short cutting point									
	70 undyed							C0069444	
1/2 circle reverse cutting needle									
	70 undyed					C0069450	C0069451		
	70 undyed				C0069460				
	70 undyed					C0069461			
Straight reverse cutting needle									
	45 undyed				C0069465				
1/2 circle round bodied needle									
	90 undyed				C0069466	C0069467			
1/2 circle round bodied needle with short cutting point									
	70 undyed							C0069475	
3/8 circle reverse cutting needle									
	70 undyed			C0069513					
3/8 circle reverse cutting needle with precision point									
	45 undyed	C0069526	C0069525	C0069524					

Novosyn®

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
3/8 circle reverse cutting needle with precision point									
 DSMP 19	70 undyed		C0069528	C0069529					
	 70 undyed			C0069529N1					
	70 undyed				C0069530				
	 70 undyed				C0069530N1				
	70 undyed					C0069531			
									
3/8 circle reverse cutting needle									
 DS 24	45 undyed			C0069534	C0069535				
									
1/2 circle round bodied needle									
 HR 26	90 undyed				C0069541	C0069540	C0069543		
									
1/2 circle round bodied needle with blunt point									
 HRN 35	70 90 undyed				C0069550	C0069549	C0069553		
									
3/8 circle reverse cutting needle with precision point									
 DSMP 24	70 undyed					C0069554			
	 70 undyed					C0069554N1			
	70 undyed			C0069555	C0069556				
	 70 undyed				C0069556N1				
	70 undyed								
									
1/2 circle round bodied needle									
 HR 40s	90 undyed						C0069559	C0069557	
									





Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle										
	90 undyed				C0069565	C0069564	C0069561	C0069562		
3/8 circle reverse cutting needle with precision point										
	90 undyed				C0069572					
1/2 circle round bodied needle										
	90 undyed					C0069585				
1/2 circle reverse cutting needle with precision point										
	45 undyed		C0069591							
1/2 circle round bodied needle										
	90 undyed				C0069594	C0069595				
	90 undyed					C0069595N1				
	90 undyed						C0069596	C0069597		
3/8 circle reverse cutting needle with precision point										
	45 undyed				C0069599					
Progressive curved reverse cutting needle with precision point										
	45 undyed			C0069600						
3/8 circle reverse cutting needle with precision point										
	45 undyed	C0069601								

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Novosyn®



ABSORBABLE SUTURES


Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
Progressive curved reverse cutting needle with precision point									
 PSMP 16 	45 undyed		C0069602						
1/2 circle reverse cutting needle with precision point									
 HSMP 12 	45 undyed	C0069606							
3/8 circle reverse cutting needle with precision point									
 DSMP 11 	45 undyed	C0069607	C0069608						
 DSMP 13f 	45 undyed	C0069617							
 DSMP 13 	45 undyed	C0069618	C0069619	C0069620					
5/8 circle reverse cutting needle									
 FS 25 	70 undyed				C0069622				
3/8 circle reverse cutting needle with precision point									
 DSMP 19 	45 undyed		C0069628	C0069629					
	45 undyed			C0069629N1					
	45 undyed				C0069630				
	45 undyed				C0069630N1				
1/2 circle reverse cutting needle with precision point									
 HSMP 15 	45 undyed			C0069636					
3/8 circle reverse cutting needle with precision point									
 DSMP 24 	45 undyed			C0069655	C0069656	C0069657			
 DSMP 39 	70 undyed				C0069659	C0069660			

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle with trocar point									
 HRT 37 	90 undyed				C0069661	C0069662	C0069663	C0069664	
 HRT 40s 	70 undyed					C0069672	C0069671	C0069670	
 HRT 40 	90 undyed						C0069673		
3/8 circle reverse cutting needle with precision point									
 2xDSMP 7 	45 undyed	C0069709							
 2xDSMP 13 	45 undyed		C0069715	C0069716					
1/2 circle round bodied needle									
 HR 22 	90 undyed				C0069720				
3/8 circle round bodied needle									
 2xDR 12 	45 undyed	C0069725							
5/8 circle round bodied needle with trocar point									
 FRT 17 	70 undyed			C0069755					
1/2 circle round bodied needle with protection point									
 HRB 48 	90 undyed							C0069772	

Novosyn®

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle with short cutting point										
	90 undyed					C0069778				
1/2 circle round bodied needle										
	45 undyed	C0069802								
1/2 circle round bodied needle with protection point										
	90 undyed						C0069822			
1/2 circle round bodied needle with short cutting point										
	90 undyed				C0069864	C0069865		C0069867		
	90 70 undyed					C0069920	C0069903			
1/2 circle round bodied needle										
	70 undyed				C0069934					
1/2 circle round bodied needle with short cutting point										
	70 undyed							C0069958		
	4 x 45 violet take off					C0088025				
	4 x 45 violet take off					C0088036				

Needle		Thread	USP (metric)							
length mm		length cm colour	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle with short cutting point										
		4 x 45 violet					C0088285			
1/2 circle round bodied needle										
		45 violet				C0088320	C0088321	C0088322		
		45 violet				C0088325				
		45 violet				C0088335				
		45 violet			C0088340					
	take off	4 x 70 violet			C0088613	C0088614				
	take off	4 x 70 violet			C0088629	C0088630				
		4 x 70 violet			C0088640	C0088641	C0088642			
		4 x 70 violet				C0088647				
	take off	4 x 70 violet				C0088660				
	take off	4 x 70 violet				C0088669				
	take off	4 x 70 violet						C0088696		

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Novosyn®

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle										
 HR 17 Easyslide take off	4 x 45 violet			C0088713	C0088714					
 HR 30ss Easyslide	70 violet						C0088716			
 HR 13 Easyslide take off	4 x 45 violet		C0088718	C0088719						
 HR 22 Easyslide take off	4 x 45 violet			C0088729	C0088730	C0088731				
 HR 22b Easyblack take off	4 x 45 violet				C0088732					
 HR 26 Easyslide take off	4 x 45 violet			C0088740	C0088741	C0088742				
 HR 30 Easyslide take off	4 x 45 violet					C0088761				
 HR 37 Easyslide take off	4 x 45 violet				C0088769	C0088770				
 HR 26s Easyslide take off	4 x 45 violet					C0088790				
 HR 37s Easyslide take off	4 x 45 violet				C0088794	C0088795				
 HR 26ss Easyslide take off	4 x 45 violet					C0088810				

Needle length mm	Thread length cm colour	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle reverse cutting needle									
 HS 26s  take off	4 x 45 violet					C0088895			
3/8 circle reverse cutting needle									
 DS 24  take off	4 x 45 violet				C0088953				
1/2 circle round bodied needle									
 HR 37  take off	4 x 45 undyed					C0089001			
 HR 17  take off	4 x 45 undyed			C0089713					
 HR 26  take off	4 x 45 undyed				C0089741				
 HR 37s  take off	4 x 45 undyed					C0089795			
1/2 circle reverse cutting needle									
 HS 26s  take off	4 x 45 undyed				C0089894				

Novosyn®

ABSORBABLE SUTURES

			Thread	USP (metric)							
			length cm colour	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)
Ligature reels											
1 thread of 250 cm			1 x 250 violet		G0058715	G0058716	G0058717	G0058719	G0058720		
1 thread of 250 cm			1 x 250 undyed		G0059715	G0059716	G0059717				
Needle			Thread	USP (metric)							
angle	length mm	diameter microns	length cm colour	8/0 (0.4)	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)
3/8 circle round bodied needle micro needle											
DRm 6			20 35 violet	G0068145 G0068146							
135°		140µ									
1/4 circle lancet needle with micro-point											
VLM 8			70 violet			G0068299					
											
3/8 circle lancet micro needle											
DLm 6			15 violet	G0068615							
140°		150µ									
1/4 circle lancet needle with micro-point											
2xVLM 6			30 violet			G0068616					
											
3/8 circle lancet micro needle											
DLm 6s			30 violet	G0068620	G0068621						
137°		200µ									
2xDLm 5			30 violet	G0068699							
140°		150µ									
2xDLm 6			45 30 violet	G0068701 G0068702							
140°		150µ									
2xDLm 6s			30 45 30 violet	G0068710	G0068711 G0068712						
137°		200µ									

Needle		Thread	USP (metric)									
angle	length mm	diameter microns	length cm colour	8/0 (0.4)	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	
3/8 circle lancet micro needle												
137°	2xDLm 6s	200μ	45 violet	G0068713								
												
1/2 circle lancet micro needle												
	2xHLm 6	150μ	30 violet	G0068715								
												
1/4 circle lancet needle with micro-point												
	VLM 10		45 violet			G0068718						
												
3/8 circle lancet micro needle												
140°	DLm 6	150μ	30 violet	G0068724								
												
1/4 circle lancet needle with micro-point												
	2xVLM 8		20 45 violet			G0068733 G0068734	G0068735	G0068736				
												
	2xVLM 6		45 violet			G0068737						
												
	2xVLM 10		45 violet			G0068738	G0068739					
												
1/2 circle lancet micro needle												
	2xHLm 8	200μ	30 45 violet	G0068740	G0068741							
												
1/4 circle lancet needle with micro-point												
	2xVLM 8		20 violet			G0068748						
												
	VLM 8		45 violet			G0068778						
												
3/8 circle round bodied needle micro needle												
135°	DRm 8	140μ	30 violet	G0068796								
												

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Novosyn®


ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)								
		8/0 (0.4)	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	
3/8 circle reverse cutting needle with precision point										
2xDSMP 6 	45 violet		G0068850							
1/2 circle lancet micro needle										
2xHLM 8 	45 violet			G0068966	G0068967	G0068968				
1/4 circle lancet needle with micro-point										
VLM 6 	45 undyed				G0069201					
1/2 circle reverse cutting needle with precision point										
HSMP 9 	45 undyed					G0069593				
HSMP 12 	45 undyed				G0069632	G0069633				
HSMP 15 	45 undyed				G0069635					
1/4 circle lancet needle with micro-point										
2xVLM 8 	45 undyed			G0069734	G0069735					
1/2 circle lancet micro needle										
2xHLM 8 	45 undyed			G0069740						
1/2 circle round bodied needle										
HR 17b  take off	8 x 45 violet							M0068153		
HR 22b  take off	8 x 45 violet					M0068182	M0068188			
HR 26b  take off	8 x 45 8 x 70 violet						M0068200 M0068201			

Needle length mm	Thread length cm colour	USP (metric)							
		8/0 (0.4)	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)
1/2 circle round bodied needle									
 HR 26s Easyslide	8 x 70 violet								M0068295
Needle length mm	Thread length cm colour	USP (metric)							
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)
 HR 26s Easyslide	8 x 70 violet						M0068305		
5/8 circle round bodied needle									
 FR 36 Easyslide	8 x 90 violet					M0088001			
1/2 circle round bodied needle									
 HR 17b Easyblack take off	8 x 45 violet		M0088011						
1/2 circle round bodied needle with short cutting point									
 HRC 26 Easyslide take off	8 x 45 violet				M0088037				
1/2 circle reverse cutting needle									
 HS 26s Easyslide take off	8 x 45 violet				M0088081				
 HS 23 Easyslide take off	8 x 45 violet			M0088090					
1/2 circle round bodied needle									
 HR 48 Easyslide take off	8 x 45 violet				M0088101	M0088102	M0088103	M0088104	

Novosyn[®]

ABSORBABLE SUTURES


Needle length mm	Thread length cm colour	USP (metric)							
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)
1/2 circle round bodied needle									
	8 x 45 violet						M0088110	M0088111	
	8 x 45 violet				M0088115	M0088116	M0088117	M0088118	
	8 x 45 violet					M0088131	M0088132		
	8 x 45 violet			M0088135	M0088136	M0088137	M0088139	M0088140	
	8 x 45 violet			M0088145	M0088146				
	8 x 45 violet				M0088160				
	8 x 45 violet				M0088164	M0088165	M0088166		
	8 x 45 violet		M0088170	M0088171	M0088172				
	8 x 45 violet		M0088177	M0088178	M0088180				
	8 x 45 violet		M0088185	M0088187					
	8 x 45 violet	M0088190	M0088191						

Needle length mm	Thread length cm colour	USP (metric)							
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)
3/8 circle reverse cutting needle									
DS 24  Easyslide take off	8 x 45 violet			M0088192					
1/2 circle round bodied needle									
HR 48  Easyslide take off	8 x 70 violet							M0088204	
HR 40s  Easyslide take off	8 x 70 violet						M0088207	M0088208	
HR 37s  Easyslide take off	8 x 70 violet				M0088212	M0088213	M0088214	M0088215	
HR 37  Easyslide take off	8 x 70 violet			M0088218	M0088219				
HR 30ss  Easyslide take off	8 x 70 violet						M0088225		
HR 26s  Easyslide take off	8 x 70 violet					M0088232	M0088233	M0088234	
HR 26  Easyslide take off	8 x 70 violet		M0088237	M0088239					
	8 x 45 undyed				M0089117				
HR 37s  Easyslide take off	8 x 45 undyed				M0089131	M0089132	M0089133		

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Novosyn®

ABSORBABLE SUTURES

Needle length mm	Thread length cm colour	USP (metric)							
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)
1/2 circle round bodied needle									
 HR 48 Easyslide take off	8 x 45 undyed					M0089140			
1/2 circle reverse cutting needle									
 HS 26s Easyslide take off	8 x 45 undyed				M0089176				
1/2 circle round bodied needle									
 HR 37s Easyslide take off	8 x 70 undyed					M0089504			
1/2 circle round bodied needle with protection point									
 HRB 48 Easyslide take off	8 x 45 undyed					M0089770	M0089771		
 HRB 37s Easyslide take off	8 x 45 undyed				M0089819	M0089820	M0089821		
 HRB 26 Easyslide take off	8 x 45 undyed			M0089828					

Handwriting practice lines consisting of 20 horizontal dotted lines.

B. Braun Surgical, S.A. | Carretera de Terrassa, 121 | 08191 Rubí | Spain
Tel. +34 935 86 62 00 | www.bbraun.com

Information according to Council Directive 93/42/EEC respectively:



Manufacturer

B. Braun Surgical, S.A. | Carretera de Terrassa, 121 | 08191 Rubí | Spain

The product trademark "Novosyn" is a registered trademark of B. Braun Melsungen AG or its affiliates.

CE 0123

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form, partial or otherwise, is not permitted.